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**STERILE and DISPOSABLE**  
**for single use only**

**Medical Device Directive (Europe)**

This device complies with the requirements of Directive 93/42 EEC concerning medical devices

Sterilized by radiation

**1. Use**

This pencil and the electrodes are intended for monopolar cutting and coagulation with various HF generator types e.g. BOWA, Erbe, Martin, Valleylab / Conmed etc.

In some cases an adapter to connect with the HF-unit will be necessary.  
The maximum permissible high frequency voltage is 6.000 V<sub>p</sub>.

**2. Operating**

- Approved use by specialist medical staff only. Use aseptic techniques when removing the product from the package.
- Remove tip protection from electrode. Check electrode for secure and tight fit in the pencil. Be also aware of this when changing the electrode.
- For 801-030, 802-030 and 802-032:  
Remove sleeve from handle and cable.
- Plug the 3-pin connector into the corresponding monopolar socket of the HF unit or use the corresponding adapter.
- For 801-030, 802-030 and 802-032:  
To activate the CUT mode, press the yellow button or switch.  
To activate the COAG mode, press the blue button or switch.  
For 803-030:  
Is activated by means of a footswitch connected with the HF unit.
- The design and construction of these products does not permit its reuse.
- BOWA ensures product compatibility only for above mentioned single-use handles and associated disposable electrodes.

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- h) Always refer to the instruction manual of the HF unit and to its general information about electrosurgical procedure.
- i) For 952-001  
Plug the handle tip through the distal opening of the suction clip until it stops and fix it in the bottom of the clip. The connecting cable of the handle can be attached to the mounting points along the tube.
- j) For 952-001 and 802-032  
Connect the end of the hose with the smoke evacuator.

### 3. Safety Warnings

- a) Sterility is only guaranteed unless package is damaged or seal is broken. Do not resterilize!
- b) Do not use the product after its expiration date.
- c) Thoroughly examine handle, cable and plug for irregularities. Do not use defective handles or electrodes.
- d) When plugging or unplugging the HF cable, always grasp the connector directly.
- e) Do not use the product in the presence of inflammable or explosive materials.

**Endogenous burns** are burns caused by high current density in the patient's tissue. Possible reasons are: the patient is inadvertently positioned in contact with electrically conductive parts. Direct contact between skin areas and HF cables may lead to electrical capacitance which in turn may cause burn.

**Exogenous burns** are those caused by the heat of ignited fluids or gases. They may also be caused by explosions. Possible reasons are among others: Ignition of skin cleaning agents and disinfectants, of narcotic gases etc.

- f) Place the active accessories away from the patient's skin when not in use. Never lay HF cables in loops.
- g) Always refer to the instruction manual of the HF unit and to its general information about electrosurgical procedures.
- h) Before operating check if the signal transmission of the monitors works without interference (e.g. without noise); this check is done by pressing the yellow button (CUT) and the blue one (COAG).
- i) **Cardiac pacemakers** may be damaged by electrosurgical current. Consult a cardiologist for the intervention. Never use electrosurgical current on outpatients fitted with pacemakers.

### 4. Environmental conditions

The allowable storage time of the original unopened packaging is 3 years after production. Make sure that the temperature range of 0°C to 40°C and a relative humidity of 0 to 75% non-condensing must be observed. During storage, exposure to sources of heat, light, sunlight and moisture is not allowed.

### 5. Disposal

Medical devices, packaging material and accessories must be disposed of according to the applicable domestic laws and regulations.