



Arthroscopy electrodes

Reusable electrodes for arthroscopy

Different use with handle or cable.

Intended use: Arthroscopy electrodes are used for electrosurgical cutting and coagulation in orthopaedics and trauma surgery.

Packaging unit: 1 pc.

Electric strength: 4250 Vp

Preparation: 75 cycles

Article no.	Article designation	Permissible combinations	Ø Shaft	Electrode type	Electrode shape
510-108	Needle electrode, arthroscopy, 90 ° angled, 1.5 x 0.8 mm	101-060, 280-035, 370-050	4 mm	Needle electrode	Angled
510-109	Needle electrode, arthroscopy, 90 ° angled, 4 x 0.8 mm	101-060, 280-035, 370-050	4 mm	Needle electrode	Angled
510-110	Knife electrode, arthroscopy, 45 ° angled, 3 x 1.5 mm	101-060, 280-035, 370-050	4 mm	Knife electrode	Cranked
510-112	Keying electrode, arthroscopy, 90 ° angled, 3 x 1.5 mm	101-060, 280-035, 370-050	4 mm	Keying electrode	Angled

Print »

BOWA-electronic GmbH & Co. KG | Heinrich-Hertz-Strasse 4-10 | 72810 Gomaringen | Germany | Tel. Fon: +49 (0) 7072-6002-0