




ClearView™


Endoscopic Ultrasound Aspiration Needle




Twist lock technology (TLT) and single hand operation for comfortable control with clear visibility.



GORE VIABIL®
BILIARY ENDOPROSTHESIS



FLEXXUS®
ENDOSCOPIC BILIARY STENT



BEAMER
ESU



To learn more about ClearView™ and other innovative products, call 800-448-6506 or visit ConMed.com.

Simplifying a Complex Procedure

CLEARVIEW™ EUS ASPIRATION NEEDLE

Designed using feedback from gastroenterologists, ClearView allows the endoscopist to perform fine needle aspiration accurately and safely for better patient outcomes.

Control at Your Fingertips

ClearView's Twist Lock Technology is crafted to allow single handed operation of the needle and sheath locks by the physician.

No more loose screws. TLT locks the pistons in place minimizing the risk of accidental needle and sheath deployment.

Visibility Through Technology

Laser etched a full 2cm from the tip of the needle.

The Power to Access Difficult Lesions

More ergonomic handle offers improved control and stability.

Variety of different puncture techniques are available to the physician.

Plastic sheath allows for increased flexibility and for obtaining aspirate in tortuous positions.



Advancing FNA

Twist Lock Technology (TLT)

- Allows physician to make one quarter turn to lock or unlock needle and sheath pistons
- Single hand operation
- Needle length adjustment (0 - 8.5 cm)
- Sheath adjustment (0-4 cm)

Ergonomically Shaped Sturdy Handle Design

- Efficient access of difficult lesions
- Offers control and stability

Treatment of the needle

- Laser etched needle for clear visibility with ultrasound
- Highly visible 2cm of treatment

Colored Plastic Sheath

- Improved visualization of sheath under endoscopic imaging
- Enhanced flexibility to access difficult lesions in angulated endoscope positions

Nitinol Stylet

- Optional extended beveled stylet to help with puncture and reduce contaminants
- Locking cap to help hold the position of the stylet during the procedure

Versions with Unique Sheath Stabilizer

- Enables physicians to precisely guide needles through larger working channels with the same small diameter needle

Enhanced Luer-Lock Design

- Secure connection to the endoscope



ClearView Ordering Information

Description	Needle Dimensions	Catalog Number
ClearView EUS Needle 19G Round	1.1mm x 2.1mm	CV2119R
ClearView EUS Needle 22G Extended Bevel	0.7mm x 1.8mm	CV1822EB
ClearView EUS Needle 22G Extended Bevel w/Sheath Stabilizer	0.7mm x 2.7mm	CV2722EB
ClearView EUS Needle 22G Round	0.7mm x 1.8mm	CV1822R
ClearView EUS Needle 22G Round w/Sheath Stabilizer	0.7mm x 2.7mm	CV2722R
ClearView EUS Needle 25G Round	0.5mm x 1.8mm	CV1825R
ClearView EUS Needle 25G Round w/Sheath Stabilizer	0.5mm x 2.7mm	CV2725R

ClearView™

Endoscopic Ultrasound Aspiration Needle



525 French Road
Utica, NY 13502

Local 315-797-8375
Toll Free 800-448-6506

ConMed.com
info@mail.conmed.com